 **BD Venflon™ Pro Safety**


BD Vialon™ Material

BD Luer-Lok™


22G 0.9 x 25mm

42 mL/min

6-SDGW34_rev08(09-20)



REF 393222


 Becton Dickinson Infusion Therapy
AB, SE-251 06 Helsingborg, Sweden

Made in Singapore bd.com

Сделано в Сингапуре

Виготовлено у Сінгапури

bd.com/symbols-glossary



STERILE R




UA.TR.116



РУ № ФСЗ 2008/03365
от 22.01.2015

CE 2797



1202077P01

LOT

2021-08-01


EXP

2024-07-31

EXP

(01) 003829 03932221



 **BD Venflon™ Pro Safety**


BD Vialon™ Material

BD Luer-Lok™


20g 1.1 x 32mm

67 mL/min

4-SDGW6_rev08(09-20)




REF 393224

 Becton Dickinson Infusion Therapy
AB, SE-251 06 Helsingborg, Sweden

Made in Singapore bd.com

Сделано в Сингапуре
Виготовлено у Сінгапури

bd.com/symbols-glossary



STERILE R




UATR.116



РУ № ФСЗ 2008/03365
от 22.01.2015




2797




1113490P45

LOT

2021-05-01




2024-04-30



(01) 003829 03932245



 **BD Venflon™ Pro Safety**


BD Vialon™ Material

BD Luer-Lok™


18G 1.3 x 32mm


103 mL/min

5-SDGW7_rev08(09-20)



REF 393226

 Becton Dickinson Infusion Therapy
AB, SE-251 06 Helsingborg, Sweden
Made in Singapore bd.com
Сделано в Сингапуре
Виготовлено у Сінгапури
bd.com/symbols-glossary




STERILE R




UATR.116



РУ № ФСЗ 2008/03365
от 22.01.2015



2797




(01) 003829 03932269


1077005P48

LOT


2021-04-01



2024-03-31





 **BD Venflon™ Pro Safety**


BD Vialon™ Material

BD Luer-Lok™


18g 1.3 x 45mm


103 mL/min

2-SDGW8_rev08(09-20)



REF 393227

 Becton Dickinson Infusion Therapy
AB, SE-251 06 Helsingborg, Sweden
Made in Singapore bd.com
Сделано в Сингапуре
Виготовлено у Сінгапури
bd.com/symbols-glossary



STERILE R




UATR.116




РУ № ФСЗ 2008/03365
от 22.01.2015




2797




1113892P41




2021-05-01




2024-04-30



(01) 003829 03932276



 **BD Venflon™ Pro Safety**


BD Vialon™ Material

BD Luer-Lok™


17G 1.5 x 45mm

133 mL/min



5-SDGW9_rev09(09-20)



REF 393228

 Becton Dickinson Infusion Therapy
AB, SE-251 06 Helsingborg, Sweden

Made in Singapore bd.com
Сделано в Сингапуре
Виготовлено у Сінгапури
bd.com/symbols-glossary



STERILE R




UATR.116



РУ № ФСЗ 2008/03365
от 22.01.2015




2797




1113937P41

LOT


2021-05-01




2024-04-30



(01) 003829 03932283



 **BD Venflon™ Pro Safety**


BD Vialon™ Material

BD Luer-Lok™


16G 1.8 x 45mm

236 mL/min


3-SDGW10_rev08(09-20)



REF 393229

 Becton Dickinson Infusion Therapy
AB, SE-251 06 Helsingborg, Sweden

Made in Singapore bd.com
Сделано в Сингапуре
Виготовлено у Сінгапури
bd.com/symbols-glossary



STERILE R




UATR.116



РУ № ФСЗ 2008/03365
от 22.01.2015



2797




1203126P41

LOT

2021-08-01




2024-07-31



(01) 003829 03932290



 **BD Venflon™ Pro Safety**


BD Vialon™ Material

BD Luer-Lok™


14G 2.0 x 45mm

270 mL/min

3-SDGW11_rev08(09-20)



REF 393230


 Becton Dickinson Infusion Therapy
AB, SE-251 06 Helsingborg, Sweden

Made in Singapore bd.com

Сделано в Сингапуре

Виготовлено у Сінгапури


bd.com/symbols-glossary



STERILE R




UATR.116



РУ № ФСЗ 2008/03365
от 22.01.2015






2797



1176670P41

2021-07-01

2024-06-30



(01) 003829 03932306